

इंटरनेट

मानक

Disclosure to Promote the Right To Information

Whereas the Parliament of India has set out to provide a practical regime of right to information for citizens to secure access to information under the control of public authorities, in order to promote transparency and accountability in the working of every public authority, and whereas the attached publication of the Bureau of Indian Standards is of particular interest to the public, particularly disadvantaged communities and those engaged in the pursuit of education and knowledge, the attached public safety standard is made available to promote the timely dissemination of this information in an accurate manner to the public.

“जानने का अधिकार, जीने का अधिकार”

Mazdoor Kisan Shakti Sangathan

“The Right to Information, The Right to Live”

“पुराने को छोड़ नये के तरफ”

Jawaharlal Nehru

“Step Out From the Old to the New”

IS 15269 (2003): Optics and Optical Instruments -
Microscopes - Magnification [PGD 22: Educational
Instruments and Equipment]



“ज्ञान से एक नये भारत का निर्माण”

Satyanarayan Gangaram Pitroda

“Invent a New India Using Knowledge”



“ज्ञान एक ऐसा खजाना है जो कभी चुराया नहीं जा सकता है”

Bhartrhari—Nitiśatakam

“Knowledge is such a treasure which cannot be stolen”

BLANK PAGE



भारतीय मानक

प्रकाशिकी और प्रकाशिक उपकरण — माइक्रोस्कोप — आवर्धन

Indian Standard

OPTICS AND OPTICAL INSTRUMENTS —
MICROSCOPES — MAGNIFICATION

ICS 37.020

© BIS 2003

BUREAU OF INDIAN STANDARDS
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG
NEW DELHI 110002

February 2003

Price Group 3

NATIONAL FOREWORD

This Indian Standard which is identical with ISO 8039:1997 'Optics and optical instruments—Microscopes—Magnification' issued by the International Organization for Standardization (ISO) was adopted by the Bureau of Indian Standards on the recommendation of the Optical and Mathematical Instruments Sectional Committee and approval of the Mechanical Engineering Division Council.

The text of ISO Standard has been approved as suitable for publication as Indian Standard without deviations. Certain conventions are, however, not identical to those used in Indian Standards. Attention is particularly drawn to the following:

- a) Wherever the words 'International Standard' appear referring to this standard, they should be read as 'Indian Standard'.
- b) Comma (,) has been used as a decimal marker while in Indian Standards the current practice is to use a point (.) as the decimal marker.

Indian Standard

**OPTICS AND OPTICAL INSTRUMENTS —
MICROSCOPES — MAGNIFICATION**

1 Scope

This International Standard specifies a series of values for the magnification of imaging components of light microscopes and defines a number of magnifying systems to which they apply.

2 Definitions

For the purposes of this International Standard, the following definitions apply.

**2.1
magnification**

act or process of changing the apparent dimensions of an object by optical methods

NOTE 1 The type of magnification such as visual or lateral should always be specified.

NOTE 2 The more general term "magnifying power" as a measure of the ability of an optical system to produce visual magnification and lateral magnification under specified operating conditions has been replaced in this International Standard by "magnification", due to the more established use of this term in practical work.

**2.2
visual magnification**

ratio of the tangents of the viewing angle of an object when observed through a magnifying system with the image at infinity to that of the object when observed by the naked eye at the reference viewing distance (250 mm)

NOTE - This ratio should be expressed in a numerical form with a multiplication sign, e.g. 10 ×.

**2.3
lateral magnification**

ratio of a given distance in the real image normal to the optical axis to the corresponding distance in the object

NOTE - This ratio should be expressed in a proportional form, e. g. 10:1.

2.4 Magnification of an objective

2.4.1 magnification of an objective with finite primary image distance

lateral magnification at the primary image formed at the distance from the objective specified in the design of the objective

2.4.2 magnification of an objective with infinite primary image distance in combination with the normal tube lens

lateral magnification at the real image produced by the combination of the objective and the 'normal' tube lens (i.e. the tube lens with which the objective is designed to operate)

NOTE — See 4.1.

2.5 tube factor

factor by which the lateral magnification at the primary image is changed by an intermediate lens or lens system inserted between the objective and the primary image

NOTE — Intermediate lenses can be fixed, interchangeable or associated with accessories having their own tube factors (see 4.2).

2.6 magnification of an eyepiece

visual magnification at the virtual image formed from the primary image by the eyepiece

NOTE — See 4.3.

2.7 projection factor

factor by which the total magnification of a microscope is changed when forming a real image of the object onto a detecting device, such as the photographic emulsion in a camera

NOTE — The image can be formed in different ways (see 4.4).

2.8 total magnification of a microscope used for visual observation

visual magnification at the virtual image formed by the microscope

NOTE — See 4.5.

2.9 total magnification of a microscope used to produce a real image

lateral magnification at the real image

NOTE — See 4.5.

3 Symbols for magnification of imaging components

Table 1 gives the symbols that shall be used when referring to the magnification of imaging components and combinations thereof and gives examples of methods of expression.

4 Methods of calculation

4.1 Magnification of an objective

The value of the magnification of an objective corrected for an infinite primary image distance is given by the ratio of the focal length of the normal tube lens to that of the objective, i. e.

$$M_{O_{\infty}} = f_{\text{NTL}}/f_{O_{\infty}}$$

where

$M_{O_{\infty}}$ is the magnification of the objective corrected for an infinite primary image distance;
 f_{NTL} is the focal length of the normal tube lens, in millimetres;
 $f_{O_{\infty}}$ is the focal length of the objective, in millimetres.

4.2 Tube factor

4.2.1 When intermediate lenses are used, the total tube factor is the product of individual factors of the intermediate lenses.

4.2.2 In the case of objectives corrected for infinite primary image distance, the value of the tube factor of a tube lens used instead of the normal tube lens, with which the objective is designed to operate, is given by the ratio of its focal length to that of the normal tube lens, i. e.

$$q = f_{\text{TL}}/f_{\text{NTL}}$$

where

q is the total tube factor;
 f_{TL} is the focal length of the tube lens, in millimetres;
 f_{NTL} is the focal length of the normal tube lens, in millimetres.

4.3 Magnification of an eyepiece

The value of the magnification of an eyepiece is the ratio of the reference viewing distance to the focal length of the eyepiece, i. e.

$$M_E = 250/f_E$$

where

M_E is the visual magnification of the eyepiece;
 f_E is the focal length of the eyepiece, in millimetres;
 250 is the reference viewing distance, in millimetres.

4.4 Projection factor

The calculation of the projection factor depends on how the image is formed.

4.4.1 If the real image is formed using a normal eyepiece intended for visual observation together with an infinitely corrected camera or projection lens focused at infinity, the value of the projection factor is given by the ratio of the focal length of the camera/projection lens to the reference viewing distance, i. e.

$$p = f_{\text{PROJ}}/250$$

where

- p is the projection factor;
- f_{PROJ} is the focal length of the camera/projection lens, in millimetres;
- 250 is the reference viewing distance, in millimetres.

4.4.2 If the real image is formed using only a normal eyepiece intended for visual observation, the value of the projection factor is given by the ratio of the distance a , from the back focal point of the eyepiece to the projected image, to the reference viewing distance, i. e.

$$p = a/250$$

where

- p is the projection factor;
- a is the distance from the back focal point of the eyepiece to the projected image, in millimetres;
- 250 is the reference viewing distance, in millimetres.

4.4.3 If the real image is formed using a lens specially designed for projection as in photomicrography, this lens can be assigned a magnification for producing the image in a given plane. This value of the magnification of the photographic projection lens, M_{PHOT} , is used to calculate the total magnification of the microscope used to produce the real image and a projection factor p is not usually applicable since the normal eyepiece is not involved in forming the image.

4.5 Total magnification

4.5.1 The value of the visual magnification of a microscope is given by the product of the magnification of the objective, the total tube factor and the visual magnification of the eyepiece, i.e.

$$M_{\text{TOT VIS}} = M_{\text{O}} \cdot q \cdot M_{\text{E}}$$

where

- $M_{\text{TOT VIS}}$ is the total (visual) magnification of the microscope;
- M_{O} is the magnification of the objective;
- q is the total tube factor;
- M_{E} is the (visual) magnification of the eyepiece.

4.5.2 The value of the total magnification of a microscope used to produce a real image using a normal eyepiece intended for visual observation, or a photographic projection lens whose projection factor has been calculated, is given by the product of the magnification of the objective, the total tube factor, the magnification of the eyepiece and the projection factor, i. e.

$$M_{\text{TOT PROJ}} = M_{\text{O}} \cdot q \cdot M_{\text{E}} \cdot p$$

where

- $M_{\text{TOT PROJ}}$ is the total (lateral) magnification of the microscope;
- M_{O} is the magnification of the objective;
- q is the total tube factor;
- M_{E} is the (visual) magnification of the eyepiece;
- p is the projection factor.

4.5.3 The value of the total magnification of a microscope used to produce a real image using a specially designed photographic lens is given by the product of the magnification of the objective, the total tube factor and the magnification of the photographic projection lens, i. e.

$M_{TOT\ PROJ} = M_O \cdot q \cdot M_{PHOT}$

where

$M_{TOT\ PROJ}$ is the total (lateral) magnification of the microscope;
 M_O is the magnification of the objective;
 q is the total tube factor;
 M_{PHOT} is the (lateral) magnification of the photographic projection lens.

5 Values and tolerances for magnification

5.1 Values for magnification

The values for the magnification of imaging components or systems shall be one of the values given in table 2. The products of quotients of any two values in the table are also to be considered as values within the table. The table may be extended by a factor of 10 per row.

Formulae for the calculation of values for specific components are given in clause 4.

5.2 Tolerances on values of magnification for specific components

Tolerances on the values of magnification shall be as given in table 3.

Table 1 — Symbols for magnification and methods of expression

System/component	Symbol	Methods of expression	
		Preferred	Alternative
Objective:			
a) corrected for finite primary image distance	M_O	$M_O = 25 : 1$	25 : 1 or 25
b) corrected for infinite primary image distance	$M_{O_{\infty}}$	$M_{O_{\infty}} = 25 \times$	25 $\times$
Eyepiece	M_E	$M_E = 10 \times$	10 $\times$
Tube factor	q	$q = 1,25 \times$	1,25
Projection factor for real image	p	$p = 0,32 \times$	0,32
Photographic projection lens	M_{PHOT}	$M_{PHOT} = 2,5 \times$	2,5 $\times$
Total magnification of microscope			
a) for visual observation	$M_{TOT\ VIS}$	$M_{TOT\ VIS} = 500 \times$	500 $\times$
b) for real image	$M_{TOT\ PROJ}$	$M_{TOT\ PROJ} = 500 : 1$	500 : 1

Table 2 — Values for magnification

				...0,32		0,4	0,5	0,63	0,8
1	1,25	1,6	2	2,5	3,2	4	5	6,3	8
10	12,5	16	20	25	32	40	50	63	80
100	125	160	200	250	320	400	500	630	800
1000	1250	1600	2000...						
NOTE 1 The values have been taken from the R10 series in ISO 3:1973, <i>Preferred numbers — Series of preferred numbers</i>. NOTE 2 The value 0,32 has been rounded from its R10 series value. NOTE 3 Besides the values in this table the following values are also in use: 1,5 - 15 - 30 - 60 - 150									

Table 3 — Tolerances on magnification

System/component	Tolerance, %
Objective	± 5
Tube factor	± 2
Projection factor	± 2
Eyepiece	± 5

Bureau of Indian Standards

BIS is a statutory institution established under the *Bureau of Indian Standards Act*, 1986 to promote harmonious development of the activities of standardization, marking and quality certification of goods and attending to connected matters in the country.

Copyright

BIS has the copyright of all its publications. No part of these publications may be reproduced in any form without the prior permission in writing of BIS. This does not preclude the free use, in the course of implementing the standard, of necessary details, such as symbols and sizes, type or grade designations. Enquiries relating to copyright be addressed to the Director (Publication), BIS.

Review of Indian Standards

Amendments are issued to standards as the need arises on the basis of comments. Standards are also reviewed periodically; a standard along with amendments is reaffirmed when such review indicates that no changes are needed; if the review indicates that changes are needed, it is taken up for revision. Users of Indian Standards should ascertain that they are in possession of the latest amendments or edition by referring to the latest issue of 'BIS Catalogue' and 'Standards: Monthly Additions'.

This Indian Standard has been developed from Doc: No. ME 31 (631).

Amendments Issued Since Publication

Amend No.	Date of Issue	Text Affected

BUREAU OF INDIAN STANDARDS

Headquarters:

Manak Bhavan, 9 Bahadur Shah Zafar Marg, New Delhi 110002
Telephones: 323 01 31, 323 3375, 323 94 02

Telegrams: Manaksanstha
(Common to all offices)

Regional Offices:

	Telephone
Central : Manak Bhavan, 9 Bahadur Shah Zafar Marg NEW DELHI 110002	323 76 17, 323 38 41
Eastern : 1/14 C.I.T. Scheme VII M, V.I.P. Road, Kankurgachi KOLKATA 700054	{ 337 84 99, 337 85 61 337 86 26, 337 91 20
Northern : SCO 335-336, Sector 34-A, CHANDIGARH 160022	{ 60 38 43 60 20 25
Southern : C.I.T. Campus, IV Cross Road, CHENNAI 600113	{ 254 12 16, 254 14 42 254 25 19, 254 13 15
Western : Manakalaya, E9 MIDC, Marol, Andheri (East) MUMBAI 400093	{ 832 92 95, 832 78 58 832 78 91, 832 78 92
Branches : AHMEDABAD. BANGALORE. BHOPAL. BHUBANESHWAR. COIMBATORE. FARIDABAD. GHAZIABAD. GUWAHATI. HYDERABAD. JAIPUR. KANPUR. LUCKNOW. NAGPUR. NALAGARH. PATNA. PUNE. RAJKOT. THIRUVANANTHAPURAM. VISAKHAPATNAM.	